


- - A soft and controlled filtration.
  - Control of flow rates , temperatures and pressures.
  - Clogging control.
  - Linear proportional dosing (peristaltic pump) of gum arabic controlled by IFS BRC validated flow meter.
  - Optimized regeneration rates.
  - Validation and traceability of sterilization cycles.
  - 3A and EHEDG guaranteed clamp mount sensors (quick to remove for inspection and calibration).
  - A communication on network: report on production , washing , history of the wine flows by carter , water flow.
  - Display of instantaneous flows.



Fauchier 4 ZAE 4 Bellevue - 33710 PUGNAC Tél: 05 57 43 65 88 fauchier@pemflow.com www.pemflow.fr/fauchier

